

# MISCELLANEOUS INSTRUMENTS



## SECTION Z



### CONTENTS

Lead Hands	1
Forceps Holders	2
Finger Ring Cutter	2
Pin Wheel	2

Note: Oxidation of lead surfaces is a normal process.



**LEAD HAND**

OM 49-0092  
9 1/2" (241mm)  
child



**LEAD HAND**

OM 49-0090  
14" (356mm)  
adult


*To prolong the life of lead hands, autoclave temperatures should not exceed 248°F (120°C). Avoid using chemical sterilization as they can cause severe corrosion.*

*Thinner and more flexible model. Longer fingers wrap on ends eliminating side tabs.*




**LEAD HAND**

OM 49-0093  
9 1/2" (241mm)  
adult  
2mm thick



**ORTHOMED FORCEPS HOLDER**

OM 39-0672    OM 39-0674    OM 39-0673  
6" (152mm)    12" (305mm)    8" (204mm)



**BUNT FORCEPS HOLDER**


OM 39-0670



**FINGER RING CUTTER**

OM 29-0430  
1/2" (165mm) chrome  
with saw blade

*Saw blade only.*



**WARTENBERG NEUROLOGICAL PIN WHEEL**

OM 28-0514 7" (178mm)