

RONGEURS



SECTION N

CONTENTS

Rongeurs



1-11



To grasp, excise tissue, degenerated disc material, and trim or sculpt bone.

ORTHO+MED
The Surgical Instrument Specialists

For soft tissue and cancellous bone only.



OM 47-0240



OM 47-0242



OM 47-0244

OM 47-0240
OM 47-0242
OM 47-0244

straight 1.5x5mm bite
30° angled 1.5x5mm bite
90° angled 1.5x5mm bite


MICRO FRIEDMAN RONGEUR

Extremely delicate
5 1/2" (152mm)

For soft tissue and cancellous bone only.




OM 47-0252
slight curve
2x5mm bite



OM 47-0253
strong curve
2x5mm bite

SYNOVECTOMY RONGEUR


Extremely delicate
6 3/4" (172mm)



OM 47-0262
5 1/2" (140mm)
curved delicate
4x11mm bite

FRIEDMAN RONGEUR

Curved delicate
5 1/2" (140mm)



OM 47-0266
30° curved
4x11mm bite


OM 47-0268
90° curved
4x11mm bite

DEAN BLUMENTHAL RONGEUR
6" (152mm)



OM 47-0267
45° curved
4x13mm bite

CHICHIRELLI RONGEUR (BLUMENTHAL)
6" (152mm)



OM 47-0263
straight delicate
2x8mm bite

OM 47-0264
curved delicate
2x8mm bite

MINI FRIEDMAN RONGEUR
6" (152mm)



OM 47-0286
curved
3x10mm bite

DEAN RONGEUR
5 3/4" (146mm)




OM 47-0274
4x9mm bite

RONGEUR #4
5 1/2" (140mm) curved



OM 47-0275
4x7mm bite


RONGEUR #4A
6 1/2" (165mm) curved



OM 47-0258
straight
3x7mm bite

OM 47-0260
curved
3x7mm bite

LEMPERT RONGEUR
6" (152mm) delicate



OM 47-0254
straight
3x7mm bite


OM 47-0256
curved
3x7mm bite

LEMPERT RONGEUR
7 1/2" (191mm) delicate



OM 47-0553
curved on side
2x10mm bite

**KLEINERT-KUTZ
SYNOVECTOMY RONGEUR**
6" (152mm) delicate




OM 47-0560

OM 47-0562

OM 47-0560
slight curve
2x13mm bite

OM 47-0562
slight curve
3x16mm bite

**KLEINERT-KUTZ
SYNOVECTOMY RONGEUR**
6" (152mm) delicate



OM 47-0558
full curve
3x8mm bite

**KLEINERT-KUTZ
SYNOVECTOMY RONGEUR**
6" (152mm) delicate



OM 47-0523
delicate

OM 47-0533

OM 47-0532

OM 47-0534

Straight
OM 47-0523 2x14 bite
OM 47-0532 3x15 bite

Curved
OM 47-0533 2x15 bite
OM 47-0534 3x15 bite

RUSKIN RONGEUR
6" (152mm)



Straight
 OM 47-0526
 4x16mm bite

OM 47-0528
 6x16mm bite

Curved
 OM 47-0527
 4x16mm bite

OM 47-0531
 6x16mm bite

RUSKIN RONGEUR
 7 1/4" (184mm)




OM 47-0525
 5x10mm bite

RUSKIN RONGEUR
 7 1/4" (184mm) curved



OM 47-0524
 5x15mm bite


ZAUFAL-JANSEN RONGEUR
 7" (178mm) curved



OM 47-0520
 3x14mm bite

OM 47-0522
 4x15mm bite

BEYER RONGEUR
 7" (178mm) curved



OM 47-0310 7" (178mm)
curved 3x16mm bite

BANE RONGEUR



OM 47-0276 6 1/2" (165mm)
curved 4x13mm bite

MEAD RONGEUR



OM 47-0285 7 1/4" (184mm)
curved 4x16mm bite

HARTMANN RONGEUR



OM 47-0308
5x15mm bite


BACON RONGEUR
8" (203mm) curved on side




Straight
OM 47-0278
8x12mm bite

Curved
OM 47-0284
8x12mm bite

LUER RONGEUR
7" (178mm)






Straight
OM 47-0302
7x15mm bite

Curved
OM 47-0306
7x15mm bite


ADSON RONGEUR
8" (203mm)

*Teeth facilitate easier
bone excision.*



OM 47-0304
with teeth
8x15mm bite

RAT TOOTH RONGEUR
8" (203mm) straight




OM 47-0480
2x12mm bite

OM 47-0482
3x12mm bite

OM 47-0484
4x12mm bite

ECHLIN DUCKBILL RONGEUR
9" (229mm) straight, curved on side










OM 47-0473
6x15mm bite

STILLE-LUER ECHLIN RONGEUR
9" (229mm) straight, curved on side











LEKSELL RONGEUR
9" (229mm), 15° curve

- | | | |
|---|---|---|
|  |  |  |
| OM 47-0492
3x16mm bite | OM 47-0486
4x16mm bite | OM 47-0487
5x16mm bite |
|  |  |  |
| OM 47-0488
6x16mm bite | OM 47-0489
7x16mm bite | OM 47-0490
8x16mm bite |




LEKSELL RONGEUR
9" (229mm), 50° curve

- | | | |
|---|---|---|
|  |  |  |
| OM 47-0494
3x13mm bite | OM 47-0495
4x13mm bite | OM 47-0497
5x13mm bite |
|  |  |  |
| OM 47-0498
6x13mm bite | OM 47-0499
7x13mm bite | OM 47-0496
8x16mm bite |



OM 47-0472
curved on side
6x13mm bite

STILLE-LUER DUCKBILL RONGEUR
9" (229mm)



Curved on side
OM 47-0461 8 1/2"

OM 47-0461
10x16mm bite



Straight
OM 47-0458 8 1/2"
OM 47-0462 10 1/2"
10x18mm bite



Curved
OM 47-0460 8 1/2"
OM 47-0464 10 1/2"
10x18mm bite


STILLE-LUER RONGEUR




**SMITH-PETERSON
LAMINECTOMY RONGEUR**
9 1/2" (241mm)



OM 47-0474
straight
3x16mm bite



OM 47-0476
curved
3x16mm bite




OM 47-0478
fully curved
3x16mm bite



DALE THORACIC RONGEUR
14 1/2" (368mm)



OM 47-0396
slight curve
3x16mm bite



SYPERT RONGEUR
14 1/2" (368mm)




OM 47-0380
4x16mm bite



OM 47-0382
5x16mm bite



OM 47-0384
8x16mm bite



OM 47-0386
12x16mm bite



OM 55-0990 5mm X 14mm
OM 55-0991 7mm X 16mm
OM 55-0992 10mm X 16mm
8 3/4" handle

ORTHO RONGEUR with PISTOL GRIP