

# WIRE AND PIN IMPLANTS

## SECTION R

### CONTENTS

Kirschner Wire	1
Steinmann Pin	2
Guide Wires	3
Shanz Pin	3
Cerclage Wire	3
Jurgan Pin Balls and General Set	4
Wire & Pin Sterilization Container	5-8
TKA Cutting Block Fixation Pins	9

Note: An internal fixation device must never be reused.

### ORTHOMED WIRE AND PIN IMPLANTS

- Precision manufactured in Germany from certified implant stainless steel.
- Longer, gently tapered points for easier penetration.
- Smoothly hand polished, causing less trauma when removed.
- K-Wire and Steinmann pins are sold six (6) to a package.

Please inquire about the availability of any sizes and styles not shown in the catalog.

#### KIRSCHNER WIRES SMOOTH

LENGTH	DIAMETER AND PRODUCT NUMBER				
	.028" (.7mm)	.035" (.9mm)	.045" (1.1mm)	.054" (1.4mm)	.062" (1.6mm)

#### DOUBLE END TROCAR

4"	OM 71-0100	OM 71-0101	OM 71-0102	OM 71-0104	OM 71-0103
5"	OM 71-0105	OM 71-0106	OM 71-0107	OM 71-0108	OM 71-0109
6"	OM 71-0021	OM 71-0022	OM 71-0023	OM 71-0024	OM 71-0025
9"	OM 71-0110	OM 71-0111	OM 71-0112	OM 71-0114	OM 71-0113
12"	OM 71-0120	OM 71-0121	OM 71-0122	OM 71-0124	OM 71-0123

#### SINGLE TROCAR Round end

4"	OM 71-0130	OM 71-0131	OM 71-0132	OM 71-0133	OM 71-0134
5"	OM 71-0170	OM 71-0171	OM 71-0172	OM 71-0174	OM 71-0175
6"	OM 71-0178	OM 71-0180	OM 71-0182	OM 71-0183	OM 71-0185
9"	OM 71-0150	OM 71-0151	OM 71-0152	OM 71-0154	OM 71-0153
12"	OM 71-0160	OM 71-0161	OM 71-0162	OM 71-0164	OM 71-0163

#### DOUBLE END DIAMOND

4"	OM 71-0200	OM 71-0201	OM 71-0202	OM 71-0203	OM 71-0204
5"	OM 71-0205	OM 71-0206	OM 71-0207	OM 71-0208	OM 71-0209
6"	OM 71-0030	OM 71-0031	OM 71-0032	OM 71-0034	OM 71-0033
9"	OM 71-0210	OM 71-0211	OM 71-0212	OM 71-0214	OM 71-0213

#### SINGLE DIAMOND Round end

4"	OM 71-0230	OM 71-0231	OM 71-0232	OM 71-0233	OM 71-0234
9"	OM 71-0240	OM 71-0241	OM 71-0242	OM 71-0244	OM 71-0243
12"	OM 71-0250	OM 71-0251	OM 71-0252	OM 71-0254	OM 71-0253

#### SUTURE PASSER SINGLE DIAMOND With suture hole

9"	—	—	OM 71-0262	—	OM 71-0263
----	---	---	------------	---	------------

#### KIRSCHNER WIRE FULLY THREADED

LENGTH	DIAMETER AND PRODUCT NUMBER			
	.035" (.9mm)	.045" (1.1mm)	.054" (1.4mm)	.062" (1.6mm)

#### DOUBLE END TROCAR

6"	OM 71-0290	OM 71-0291	OM 71-0292	OM 71-0293
9"	OM 71-0301	OM 71-0302	—	OM 71-0303

#### SINGLE TROCAR

9"	OM 71-0304	OM 71-0305	—	OM 71-0310
----	------------	------------	---	------------

#### DOUBLE END DIAMOND

9"	OM 71-0311	OM 71-0312	—	OM 71-0313
----	------------	------------	---	------------


#### SINGLE DIAMOND

9"	OM 71-0314	OM 71-0315	—	OM 71-0307
----	------------	------------	---	------------



Bend K-Wire up to .062" in diameter easily with our Gratloch™ Wire Bender OM 48-0370. See Section S.

Drawings shown on this page are not actual size.



To aid in alignment, reduction, and fixation of fractured bone.



STEINMANN PINS SMOOTH

LENGTH	DIAMETER AND PRODUCT NUMBER							
	5/64" (2.0mm)	3/32" (2.4mm)	7/64" (2.8mm)	1/8" (3.2mm)	9/64" (3.5mm)	5/32" (4.0mm)	.177" (4.5mm)	3/16" (4.8mm)

DOUBLE END TROCAR 

9"	OM 71-00450	OM 71-0453	OM 71-0456	OM 71-0459	OM 71-0462	OM 71-0465	OM 71-0468	OM 71-0469
----	-------------	------------	------------	------------	------------	------------	------------	------------

SINGLE TROCAR 3 shank end 

9"	OM 71-00480	OM 71-0483	OM 71-0486	OM 71-0489	OM 71-0492	OM 71-0495	OM 71-0498	—
----	-------------	------------	------------	------------	------------	------------	------------	---

SINGLE TROCAR Round end 

9"	OM 71-00510	OM 71-0513	OM 71-0516	OM 71-0519	OM 71-0522	OM 71-0525	OM 71-0528	OM 71-0529
12"	OM 71-00674	OM 71-0676	OM 71-0677	OM 71-0680	OM 71-0681	OM 71-0682	OM 71-0683	—

DOUBLE END DIAMOND 

9"	OM 71-00600	OM 71-0603	OM 71-0606	OM 71-0609	OM 71-0612	OM 71-0615	OM 71-0618	OM 71-0619
----	-------------	------------	------------	------------	------------	------------	------------	------------

SINGLE DIAMOND 3 shank end 

9"	OM 71-00630	OM 71-0633	OM 71-0636	OM 71-0639	OM 71-0642	OM 71-0645	OM 71-0648	—
----	-------------	------------	------------	------------	------------	------------	------------	---

SINGLE DIAMOND Round end 


9"	OM 71-00660	OM 71-0663	OM 71-0666	OM 71-0669	OM 71-0672	OM 71-0675	OM 71-0678	OM 71-0679
----	-------------	------------	------------	------------	------------	------------	------------	------------

SINGLE DIAMOND With suture hole, 3 shank end 

9"	OM 71-00880	OM 71-0883	OM 71-0886	OM 71-0889	OM 71-0892	OM 71-0895	OM 71-0898	—
----	-------------	------------	------------	------------	------------	------------	------------	---

STEINMANN PINS FULLY THREADED

LENGTH	DIAMETER AND PRODUCT NUMBER							
	5/64" (2.0mm)	3/32" (2.4mm)	7/64" (2.8mm)	1/8" (3.2mm)	9/64" (3.5mm)	5/32" (4.0mm)	.177" (4.5mm)	3/16" (4.8mm)

DOUBLE END TROCAR 

9"	OM 71-00750	OM 71-0753	OM 71-0756	OM 71-0759	OM 71-0762	OM 71-0765	OM 71-0768	OM 71-0769
----	-------------	------------	------------	------------	------------	------------	------------	------------

SINGLE TROCAR Round end 

9"	OM 71-00770	OM 71-0773	OM 71-0776	OM 71-0779	OM 71-0782	OM 71-0785	OM 71-0788	OM 71-0789
----	-------------	------------	------------	------------	------------	------------	------------	------------

DOUBLE END DIAMOND 

9"	OM 71-00730	OM 71-0733	OM 71-0736	OM 71-0739	OM 71-0742	OM 71-0745	OM 71-0748	OM 71-0749
----	-------------	------------	------------	------------	------------	------------	------------	------------

SINGLE DIAMOND Round end 

9"	OM 71-00700	OM 71-0703	OM 71-0706	OM 71-0709	OM 71-0712	OM 71-0715	OM 71-0718	OM 71-0719
----	-------------	------------	------------	------------	------------	------------	------------	------------

Drawings shown on this page are not actual size.

STEINMANN PINS CALIBRATED (GUIDE WIRES) with 1/2" graduations

LENGTH	DIAMETER AND PRODUCT NUMBER			
	5/64" (2.0mm)	3/32" (2.4mm)	7/64" (2.8mm)	1/8" (3.2mm)

SINGLE TROCAR 3 shank end 

9"	—	OM 71-0833	—	—
----	---	------------	---	---

SINGLE DIAMOND 3 shank end 

9"	—	OM 71-0853	OM 71-0856	OM 71-0859
----	---	------------	------------	------------

SINGLE DIAMOND Round end 

9"	—	OM 71-0873	—	—
----	---	------------	---	---

SCHANZ PINS (Partially Threaded)

PRODUCT NUMBER	OVERALL LENGTH	THREAD LENGTH	DIAMETER
OM 72-0220	6"	45mm	5mm
OM 72-0222	6 3/4"	45mm	5mm
OM 72-0230	7 1/2"	45mm	5mm
OM 72-0224	8"	45mm	5mm
OM 72-0226	10"	45mm	5mm

CERCLAGE WIRE 10 meters in length

PRODUCT NUMBER	DIAMETER	APPROX. GAUGE
OM 72-0102	0.2mm	36
OM 72-0103	0.3mm	30
OM 72-0104	0.4mm	27
OM 72-0105	0.5mm	25
OM 72-0106	0.6mm	23
OM 72-0107	0.7mm	22
OM 72-0108	0.8mm	21
OM 72-0109	0.9mm	20
OM 72-0110	1.0mm	19
OM 72-0112	1.2mm	18
OM 72-0115	1.5mm	17
OM 72-0116	1.8mm	13





**JURGAN REFILL RACKS  
STANDARD WHITE**

W-Series pin balls

Six balls of one size per rack

OM JR-W028	.028" (.7mm)
OM JR-W035	.035" (.9mm)
OM JR-W045	.045" (1.1mm)
OM JR-W062	.054" -.062" (1.6mm)
OM JR-W564	5/64" (2mm, .065")
OM JR-W332	3/32" (2.5mm)

**JURGAN REFILL RACKS  
COLOR CODED**

W-Series pin balls

Six balls of one size per rack

(red, dark blue, yellow, green, blue, cream)

OM JR-W028-RD	.028" (.7mm)
OM JR-W035-DB	.035" (.9mm)
OM JR-W045-YL	.045" (1.1mm)
OM JR-W062-GN	.054" -.062" (1.6mm)
OM JR-W564-BL	5/64" (2mm, .065")
OM JR-W332-CR	3/32" (2.5mm)



**JURGAN SCREWDRIVER**

Standard - for use with W and X series balls  
OM JR-S200



**JURGAN PIN BALL GENERAL SET**

For capping K-Wires  
OM JR-C50 standard white

Set includes:  
OM JR-W045 .045" (1.1mm)  
OM JR-W062 .054" -.062" (1.6mm)  
OM JR-W564 5/64" (2mm, .065")  
OM JR-W332 3/32" (2.5mm)  
OM JR-S200 Screwdriver  
OM JR-C10 Sterilization Case

*Designed to provide easy access  
and identification.*



**ORTHOMED K-Wire AND STEINMANN  
PIN STERILIZATION CASE**

OM 71-0010  
11 3/4" x 5 1/4" x 2" (closed)  
(298mm x 133mm x 51mm)

The OrthoMed K-Wire and Steinmann Pin Sterilization Case holds popular 4" (102mm), 6" (152mm) and 9" (229mm) K-Wires and 9" (229mm) Steinmann Pins.

Accommodates four each of the following K-Wire diameters:

.028"	(.7mm)
.035"	(.9mm)
.045"	(1.1mm)
.062"	(1.6mm)


Accommodates two each of the following Steinmann Pin diameters:

5/64"	(2.0mm)
3/32"	(2.4mm)
7/64"	(2.8mm)
1/8"	(3.2mm)
9/64"	(3.5mm)
5/32"	(4.0mm)
.177"	(4.5mm)

- Removable lid fits under case – taking up less room in the O.R.
- Nylon coated brackets prevent metal-to-metal contact – protecting sharp ends.
- Contents are held in place when closed – preventing shifting.
- Secure locking side brackets help prevent accidental opening.
- Handles at each end for easy transport.
- Anodized aluminum case is lightweight and built to withstand abuse.
- Fully autoclavable to 270°F (132°C).



*Case in closed position.*



The OrthoMed, Inc. K-Wire and Steinmann Pin Sterilization Case holds popular 9" K-Wires and Steinmann Pins.

Accommodates four each of the following K-Wire diameters:

- .028" (.7mm)
- .035" (.9mm)
- .045" (1.1mm)
- .062" (1.6mm)

Accommodates four each of the following Steinmann Pin diameters:

- 5/64" (2.0mm)
- 3/32" (2.4mm)
- 7/64" (2.8mm)
- 1/8" (3.2mm)
- 9/64" (3.5mm)
- 5/32" (4.0mm)
- .177" (4.5mm)

- Removable lid fits under case – taking up less room in the O.R.
- Nylon coated brackets prevent metal-to-metal contact – protecting sharp ends.
- Contents are held in place when closed – preventing shifting.
- Secure locking side brackets help prevent accidental opening.
- Handles at each end for easy transport.
- Anodized aluminum case is lightweight and built to withstand abuse.
- Fully autoclavable to 270°F (132°C).



Case in closed position.

**ORTHOMED K-Wire AND STEINMANN PIN STERILIZATION CASE**

OM 71-0009  
13 5/8" x 5 1/4" x 2 1/4" (closed)  
(346mm x 133mm x 57mm)



Designed to provide easy access and identification.

Case in open upright position

The OrthoMed, Inc. K-Wire sterilization container holds 4" or 6" K-Wires.

Accommodates six each of the following K-Wire diameters:

- .028" (.7mm)
- .035" (.9mm)
- .045" (1.1mm)
- .054" (1.4mm)
- .062" (1.6mm)

- K-Wire sizes are clearly noted on the inside of case for quick identification and reordering.
- Lid folds back and acts as a support to hold the case in an upright position – taking up less room in the O.R.
- Top of case locks when closed helping to prevent accidental opening.
- Case has a one piece design – eliminates the possibility of losing or switching case parts.
- Anodized aluminum case is lightweight and built to withstand abuse.
- Fully autoclavable to 270°F.



Case in closed position.

**ORTHOMED K-Wire STERILIZATION CONTAINER**

OM 71-0008  
6 1/4" x 4 1/8" x 3/4" (closed)  
(158mm x 105mm x 19mm)





The OrthoMed, Inc. K-Wire and Steinmann Pin Sterilization Case holds popular 9" K-Wires and Steinmann Pins.

Accommodates twelve each of the following K-Wire diameters:

- .035" (.9mm)
- .045" (1.1mm)
- .062" (1.6mm)

Accommodates four each of the following Steinmann Pin diameters:

- 5/64" (2.0mm)
- 3/32" (2.4mm)

Accommodates two each of the following Steinmann Pin diameters:

- 7/64" (2.8mm)
- 1/8" (3.2mm)
- 9/64" (3.5mm)
- 5/32" (4.0mm)

- Removable lid fits under case – taking up less room in the O.R.
- Nylon coated brackets prevent metal-to-metal contact – protecting sharp ends.
- Contents are held in place when closed – preventing shifting.
- Secure locking side brackets help prevent accidental opening.
- Handles at each end for easy transport.
- Anodized aluminum case is lightweight and built to withstand abuse.
- Fully autoclavable to 270°F (132°C).



Case in closed position.

**ORTHOMED K-Wire AND STEINMANN PIN STERILIZATION CASE**

OM 71-0011  
13 5/8" x 5 1/4" x 2 1/4" (closed)  
(346mm x 133mm x 57mm)



**BONE PIN with 3 Flutes**

OM 72-0300 3.5" 3.2mm diameter



**BONE PIN Smooth**

OM 72-0360 3" (76mm) round end, 3.2mm diameter  
OM 72-0364 3.5" (90mm) square end, 3.2mm diameter



**TROCAR TAP**

OM 72-0330 3.2 mm Dia 45 mm Length



**TROCAR TAP**

OM 72-0346 3.2 mm DIA 60 mm Length



**TAP**

OM 72-0335 3.2 DIA 70 mm Length 2 Flutes



**TAP**

OM 72-0338 3.2mm X 67mm



**JIG PIN**

OM 72-0366 3.2mm x 76mm smooth tri end