

CARTILAGE AND TENDON INSTRUMENTS

SECTION P




To retrieve, grasp, hold, cut and separate cartilage and tendons.

CONTENTS


Cartilage Instruments	1-2
Tendon Instruments	2-3
Tendon Grasping Forceps	4
Tendon Strippers	5
Cartilage Knives	5
Tissue Graspers	6




MARTIN CARTILAGE SCISSORS



MUELLER CAPSULE SCISSORS



DINGMAN FORCEPS




MARTIN CARTILAGE CLAMP




OM 46-0818
8" (203mm)
slightly curved

WALTON CARTILAGE CLAMP




OM 46-0744
5 3/8" (137mm)
serrated, straight

MENISCUS CLAMP




OM 46-0746
5 3/8" (137mm)
serrated, angled 20°



*For grasping tags of
the meniscus.*

OM 46-0748
7 1/4" (184mm)
serrated jaw with 1x2 tip

**PRESTON LIGAMENTUM FLAVUM
FORCEPS (MENISCUS)**



OM 46-0736 6" (152mm) angular shaft
OM 46-0737 7 1/2" (191mm) angular shaft
OM 46-0738 9 1/2" (241mm) angular shaft

BRAND TENDON PULLING FORCEPS



OM 90-1883 sharp
OM 90-1884 blunt

MENISCAL CLAMP 7"



OM 46-0739 8" (203mm)
serrated
with sharp teeth
OM 46-0735 8" (203mm)
serrated
with sharp teeth
flexible shaft

KLEINERT-KUTZ TENDON RETRIEVER

INLET CLOSURE NON-DISPOSABLE

OM 1700040
9 1/2" (241mm)
with cleaning port




CAROLL TENDON PULLING FORCEPS

OM 46-0740
5" (127mm)
curved shaft
serrated



SUTURE GUIDE

OM 6060 10mm
OM 6070 12mm



HOEN GRASPING FORCEPS

*Designed with one long
jaw. Longer jaw extends
5mm and pierces the
host tendon. The jaw then
grasps the graft. When
pulled back it draws the
graft through the host to
complete the first weave.*



OM 46-0742 straight
OM 46-0743 curved

TENDON WEAVING FORCEPS

6 1/4" (159mm)



OM 46-0730
3mm

OM 46-0731
4mm


OM 46-0732
5mm

OM 46-0733
6mm

OM 46-0734
7mm

BUNNELL TENDON STRIPPER

6" (152mm)



OM 46-0766
curved left


OM 46-0770
curved right

*Also known as a
Meniscus Knife.
Useful for
fasciotomies.*

OM 46-0762
straight
chisel type


SMILLIE KNIFE

6 3/4" (172mm)



OM 46-0741
10" (254mm)
straight edge with guards
curved blade

LOWE-BRECK CARTILAGE KNIFE



OM 46-0772
10" (254mm)
concave edge
with guards curved blade

DOWNING CARTILAGE KNIFE



CEMENT GRASPER

OM RW8800 9 3/4" overall length (7" shaft, 15mm x 3.5mm cup)



INTRA-ARTICULAR TISSUE GRASPER/RONGEUR

OM 90-1790-01 8" overall length (5" shaft length)
OM 90-1790-03 10" overall length (7" shaft length)
OM 90-1790-02 12" overall length (9" shaft length)



WILLIAMS IVD RONGEUR

5 1/4" (135mm)
OM RW4849
2mm with 1X2 teeth