


FOOT AND ANKLE



SECTION D

CONTENTS

Nail Splitters	1-2
Tissue Nippers	2
Nail Nippers	3-5
Osteotomy Guide	5
Nail Forceps	5-6
Nail Scissors	6
Hand Instruments	7-9
Dermal Instruments	10-11
Calcaneal/Metatarsal Spreaders	12



To split, cut, or remove finger and toe nails.



OM 69-0320
4" (102mm)
stainless



OM 69-0322
4 1/2" (114mm)
stainless




OM 69-0326
5" (127mm)
heavy jaw
stainless



OM 69-0327
5" (127mm)
narrow jaw
stainless



OM 69-0340
6" (152mm)
stainless



NAIL SPLITTER

Note: Replacement springs are available to order. Refer to the corresponding product number.



OM 69-0341
4" (102mm)
tapered
stainless



OM 69-0342
4 1/2" (114mm)
tapered
stainless



OM 69-0343
5" (127mm)
tapered
stainless



OM 69-0344
6" (152mm)
tapered
stainless*




OM 69-0448
4 1/2" (114mm)
delicate
stainless*



NAIL SPLITTER

*Ideal for use in P/A

Note: Our splitters have the true feel of superb craftsmanship. The cutting edge is expertly ground and specially heat treated. For a sure grip, our stainless splitters come with grooved handles.



OM 69-0930
5" (127mm)




OM 69-0932
6" (152mm)

Short jaw design provides increased leverage for easier cutting of thick nails. Less strain on hand and splitter.

SHORT JAW NAIL SPLITTER

Stainless



OM 69-0290
4" (102mm)
6mm stainless



OM 69-0291
4" (102mm)
9mm stainless




OM 69-0292
4 1/2" (114mm)
11mm stainless




OM 69-0294
5" (127mm)
14mm stainless

TISSUE NIPPER



NAIL NIPPER

Heavy




OM 69-0430
5 1/2" (140mm)
concave
stainless




OM 69-0432
5 1/2" (140mm)
straight
stainless



OM 69-0438
6 1/2" (165mm)
concave
super size
stainless



Double action helps prevent hand fatigue. Ideal for mycotic and onychauxis nails.




OM 69-0445
straight



OM 69-0446
concave

MYCOTIC NAIL NIPPER

6" (152mm) stainless double action
barrel spring



OM 69-0418
5 1/2" (140mm)
concave
narrow
stainless


NAIL NIPPER



OM 69-0410
4 1/2" (114mm)
concave
chrome
barrel spring

OM 69-0412
5 1/2" (140mm)
concave
chrome
barrel spring

NAIL NIPPER



OM 69-0440
OM 69-0442



OM 69-0441

Angled concave
OM 69-0440 stainless, sharp end
OM 69-0441 stainless, blunt end
OM 69-0442 chrome, sharp end

NAIL NIPPER

5 1/2" (140mm)




OM 69-0258
4" (114mm)
6mm
chrome
ring handle

CUTICLE NIPPER



OM 69-0998
5" (125mm) arm
1.75" (45mm) blade

OSTEOTOMY GUIDE




*Very delicate anvil
design slides easily
under nails.*



OM 69-0348
5" (127mm)
anvil pattern
narrow blade
stainless

**INGROWN NAIL
FORCEPS**



OM 69-0792



OM 69-0793

OM 69-0792
3mm
narrow

OM 69-0793
6mm
broad

PLATYPUS NAIL REMOVING FORCEPS



5 1/2" (140mm)



*Used to grasp nail
spicules for removal.*

OM 69-0352
5 1/2" (140mm)


**NAIL SPICULE
FORCEPS**



Straight
OM 69-0820


Curved
OM 69-0824

MOLESKIN SCISSORS
7 1/2" (191mm) One blade serrated



OM 69-0805
6" (152mm)
one serrated blade

**INGROWN NAIL
SPLITTING SCISSORS**




Curved
OM 69-0106
3 1/2" (89mm)

NAIL SCISSORS




HAND INSTRUMENTS, DOUBLE ENDED SOLID OCTAGON HANDLES




CURETTE EXCAVATOR

OM 69-0488
5 1/2" (140mm)
double ended
1.5 and 2.0mm with hole



CURETTE EXCAVATOR

OM 69-0492
5 1/2" (140mm)
double ended
1.5 and 2.5mm with hole




PROBE AND PACKER

OM 69-0495
6" (152mm)
double ended




SPATULA AND PACKER

OM 69-0498
5 3/4" (146mm)
double ended




EXCAVATOR AND PACKER

OM 69-0502
5 3/4" (146mm)
double ended




EXCAVATOR AND PROBE

OM 69-0518
6" (152mm)
double ended




NAIL RASP

OM 69-0538
6 3/4" (172mm)
double ended, angled




JOHNSON-TUCKER HOOK AND FORK

OM 69-0790
6 1/2" (165mm)
double ended, delicate



HAND INSTRUMENTS, SINGLE ENDED SOLID OCTAGON HANDLES



CURETTE EXCAVATOR (VERUCCA)

OM 69-0457
5", 4mm without hole
single ended

*With McGlamry
Bullneck, Will resist
bending.*

1.0mm
1.5mm

2.0mm
2.5mm



CURETTE EXCAVATOR

OM 69-0457
5", 4mm without hole
single ended

OM 69-0461 5", 1.5mm
OM 69-0465 5", 2.5mm
OM 69-0463 5", 2.0mm

1.0mm
1.5mm


2.0mm
2.5mm



CURETTE EXCAVATOR


OM 69-0468
5", 1.0mm with hole
single ended

OM 69-0478 5", 1.5mm
OM 69-0484 5", 2.5mm
OM 69-0480 5", 2.0mm




NUCLEUS KNIFE

OM 69-0524
5 1/4" small
single ended




OM 69-0528
5 1/4" large
single ended




BONE FILE #92A

OM 69-0530
6" (152mm) single ended




BONE FILE #92B

OM 69-0534
6" (152mm) single ended




NAIL SPLITTER

OM 69-0557 5" small blade straight, single ended
OM 69-0562 5" small blade angled, single ended




NAIL SPLITTER

OM 69-0560 5" large blade straight, single ended
OM 69-0568 5" large blade angled, single ended




LOCKE ELEVATOR

OM 69-0579
4 1/2" (114mm) narrow blade
single ended




OM 69-0590
5" (127mm) wide blade
single ended

*With McGlamry
Bullneck, Will resist
bending.*

**GOLF STICK KNIFE #41**


OM 69-0618
5 1/2" (140mm)
single ended

**GOLF STICK KNIFE #43**


OM 69-0620
5 1/2" (140mm)
single ended

**INGROWN NAIL INSTRUMENT**


OM 69-0595
5" (127mm)
single ended

**SOFT CORN KNIFE**

OM 69-0622
5" (127mm)
round blade
straight, single ended

**SOFT CORN SPOON KNIFE**


OM 69-0623
5" (127mm) oval blade
slightly curved, single ended

**CORN KNIFE**


OM 69-0625
5 1/2" (140mm) curved
single ended

**CHISEL**

OM 69-0630
5 1/4" (133mm) medium blade
straight edge
single ended

**CHISEL**


OM 69-0640
5 1/4" (133mm) large blade
curved edge
single ended

**CHISEL**

OM 69-0645
5 1/4" (133mm) extra large blade
curved edge
single ended

**CHISEL**


OM 69-0650
5 1/4" (133mm) spade pattern
curved edge
single ended



OM 29-0580
3 3/4" (95mm)
standard
crimped loops

OM 29-0582
3 3/4" (95mm)
extra fine
crimped loops

**SCHAMBERG COMEDONE
EXTRACTOR**



OM 29-0584
3 3/4" (95mm)
square loops

**SCHAMBERG COMEDONE
EXTRACTOR**




OM 29-0592
5 3/4" (146mm)
fenestrated spoon
ends

**UNNA COMEDONE
EXTRACTOR**



OM 29-0604
6 1/4" (159mm)
one fenestrated cup
and curved lancet
with cap

**ORTHOMED COMEDONE
EXTRACTOR**



OM 29-0606
5" (127mm)
one fenestrated cup
and lancet blade with
cap

**GREG COMEDONE
EXTRACTOR**



OM 29-0566
4 1/2" (114mm)
with lancet
and one fenestrated
cup

**SAALFELD COMEDONE
EXTRACTOR**



OM 29-0511 2mm
OM 29-0512 3mm
OM 29-0513 4mm
OM 29-0514 5mm
OM 29-0515 6mm
OM 29-0517 8mm

KEYES PUNCH


4" (102mm)



OM 29-0495 1mm
OM 29-0496 2mm
OM 29-0497 3mm
OM 29-0498 4mm
OM 29-0499 5mm
OM 29-0500 6mm

FOX CURETTE

5 1/2" (140mm) round



OM 29-0503



OM 29-0505




OM 29-0507

OM 29-0503 #1 small
OM 29-0505 #2 medium
OM 29-0507 #3 large

PIFFARD CURETTE

5 1/2" (140mm) oval




OM 35-0515
4" (102mm)
angled on flat
sharp outside edges

BIRO DERMAL NAEVUS SCISSORS

Provides excellent joint exposure without blocking intra-articular or osteotomy access. Helps prevent slippage or falling out of the joint by placing the arms on either side of the area to be distracted, driving two K-Wires and opening the joint.

OM 48-0223
6 3/4" (171mm)
for up to .062"
(1.6mm)

OM 48-0224
6 3/4" (171mm)
for up to 7/64"
(2.8mm)



JOINT AND CALCANEAL SPREADER




OM 90-1880
7" (178mm)
Smooth pads 15mm x 12mm


OM 90-1881
7" (178mm)
Grooved pads 15mm x 12mm

CALCANEAL SPREADER

OM 69-0990
4 1/2" (114mm)
8mm width X 28mm
depth blades



COX METATARSAL SPREADER



Note:

- Specially designed thin serrated blades fit anatomically correct against the metatarsals.
- Precise design and German manufacture assure slip-proof hold on metatarsals.
- Affords excellent visibility for the surgeon.
- Designed by Dr. Kenneth Cox.