


GENERAL INSTRUMENTS

SECTION A



To dissect, hold, clamp, manipulate and suture soft tissue.

ORTHOMED
The Surgical Instrument Specialists

CONTENTS

Super Cut Scissors	1-3
Operating Scissors	4-7
Dissecting Scissors	4-9
Suture scissors	9-10
Tungsten Carbide Scissors	11-12
Bandage Scissors	13
Utility Scissors/Wire Scissors	14-15
Hemostatic Forceps	15-19
Intestinal Forceps	20-21
Tubing and Utility Forceps	21-20
Towel and Sponge Forceps	22-23
Dressing & Tissue Forceps	23-27
Splinter Forceps	28-29
Needle Holders	29-32
Suction Tubes	33-34
Scalpel Handles	34-35
Probes/Suture Instruments	36-38
Bi-Polar Instruments	38


SUPER-CUT SCISSORS

Super-Cut Scissors represent a truly remarkable advance in scissors technology. Using the finest surgical steel, the latest heat treating methods and master craftsmanship, the upper blade of the scissors is ground to a razor-sharp edge. This process creates scissors with an unprecedented smooth and effective cutting action.

Benefits of Super-Cut Scissors:

- Cut through soft and fibrous tissue without effort.
- Superb control permits precision in dissection.
- Less hand fatigue for the surgeon and less trauma for the patient.
- Retain sharpness longer than ordinary scissors.
- Ring handles are black-coated for easy identification.
- Fully autoclavable.
- Backed by OrthoMed's lifetime guarantee.
- Not recommended for fiberglass or plaster splints.

Surgeons who use our Super-Cut Scissors are convinced of their superior performance and won't do without them on their surgical trays.



SUPER-CUT MAYO SCISSORS

Beveled blades, blunt points


Straight

OM 35-0280 5 1/2" (140mm)
OM 35-0292 6 3/4" (171mm)
OM 35-0316 9" (229mm)

Curved

OM 35-0286 5 1/2" (140mm)
OM 35-0298 6 3/4" (171mm)
OM 35-0322 9" (229mm)
OM 35-0340 11" (279mm)

Flat, blunt scissors for dissecting soft tissue.



Straight

OM 35-0364
4 1/2" (114mm)
blunt points

Curved

OM 35-0370
4 1/2" (114mm)
blunt points



SUPER-CUT METZENBAUM SCISSORS

Blunt points

OM 35-0382

OM 35-0394

OM 35-0396

OM 35-0406

OM 35-0418

OM 35-0430

Straight

OM 35-0380 6" (152mm)
OM 35-0392 7" (171mm)

Curved

OM 35-0382 6" (152mm)
OM 35-0394 7" (178mm)
OM 35-0396 7" (178mm)
OM 35-0406 8" (203mm)
OM 35-0418 9" (229mm)
OM 35-0430 11" (279mm)


delicate



**SUPER-CUT
IRIS SCISSORS**

Sharp points

*Fine sharp/sharp
blades for cutting
delicate tissue
during plastic,
hand, and minor
vascular surgery.*



OM 35-0450

OM 35-0452

OM 35-0454

OM 35-0460

Straight
OM 35-0450 3 1/2" (89mm) delicate
OM 35-0454 4 1/2" (114mm)

Curved
OM 35-0452 3 1/2" (89mm) delicate
OM 35-0460 4 1/2" (114mm)

*Delicate short
narrow tips for
dissecting fine
tissue.*




**SUPER-CUT
STEVENS SCISSORS**

Blunt points

Straight
OM 35-0484 4 1/2" (114mm)

Curved
OM 35-0490 4 1/2" (114mm)
OM 35-0494 5 1/4" (133mm)



**SUPER-CUT
JAMISON SCISSORS**

6 1/4" (159mm)

*Very delicate
blunt scissors
with long handles
for fine tissue
dissection.*



OM 35-0496
curved blunt points
very delicate

OrthoMed developed these remarkably sharp scissors in response to the demands of doctors, nurses, and other O.R. personnel.

Hospital tests show that these scissors cut multiple layers of stockinette bandages, drapes, and moleskin with ease.

- Very little strength is required for cutting, putting less strain on thumb and forearm.
- Perfect cutting control because you don't have to force them to cut.
- Scissors cut faster, saving critical time in the O.R.
- Retain their sharpness longer than ordinary scissors.
- Not recommended for cutting plaster.

The delicate offset shanks and the sharp Super-Cut edge provide better leverage. The fine knob slides easily under finger and toe bandages.



OM 35-0202
5 1/2" (140mm)
shanks angled on side

SUPER-CUT HI-LEVEL BANDAGE SCISSORS

Cuts bandages and multiple layers of dressing with ease. Razor-sharp blades reduce hand fatigue.



OM 35-0210 OM 35-0212

OM 35-0210
5 1/2" (140mm)
stainless

OM 35-0212
7 1/4" (184mm)
stainless


SUPER-CUT LISTER BANDAGE SCISSORS

All the benefits of our Super-Cut Lister Scissors plus an oversized finger ring, which helps relieve pressure on the third metacarpal.



OM 35-0214
8" (203 mm)
stainless


SUPER-CUT LISTER BANDAGE SCISSORS WITH LARGE FINGER RING




Straight
OM 35-0220
5 1/2" (140mm)
sharp/blunt

Curved
OM 35-0250
5 1/2" (140mm)
sharp/blunt

SUPER-CUT OPERATING SCISSORS



OPERATING SCISSORS



OM 33-0470 OM 33-0472




OM 33-0474

Straight

OM 33-0300 4 1/2" (114mm) sharp/blunt
OM 33-0302 4 1/2" (114mm) sharp/sharp
OM 33-0304 4 1/2" (114mm) blunt/blunt
OM 33-0310 5" (127mm) sharp/blunt
OM 33-0312 5" (127mm) sharp/sharp
OM 33-0314 5" (127mm) blunt/blunt
OM 33-0320 5 1/2" (140mm) sharp/blunt
OM 33-0322 5 1/2" (140mm) sharp/sharp
OM 33-0324 5 1/2" (140mm) blunt/blunt
OM 33-0370 6 1/2" (165mm) sharp/blunt
OM 33-0372 6 1/2" (165mm) sharp/sharp
OM 33-0374 6 1/2" (165mm) blunt/blunt


Curved

OM 33-0470 4 1/2" (114mm) sharp/blunt
OM 33-0472 4 1/2" (114mm) sharp/sharp
OM 33-0474 4 1/2" (114mm) blunt/blunt
OM 33-0482 5" (127mm) sharp/blunt
OM 33-0484 5" (127mm) sharp/sharp
OM 33-0486 5" (127mm) blunt/blunt
OM 33-0490 5 1/2" (140mm) sharp/blunt
OM 33-0492 5 1/2" (140mm) sharp/sharp
OM 33-0494 5 1/2" (140mm) blunt/blunt
OM 33-0520 6 1/2" (165mm) sharp/blunt
OM 33-0522 6 1/2" (165mm) sharp/sharp
OM 33-0524 6 1/2" (165mm) blunt/blunt



PLASTIC SURGERY SCISSORS

4 3/4" (121mm)



OM 33-0532 OM 33-0538

Straight


OM 33-0532 sharp/sharp
OM 33-0534 sharp/blunt
OM 33-0536 blunt/blunt

Curved

OM 33-0538 sharp/sharp
OM 33-0540 sharp/blunt
OM 33-0542 blunt/blunt

Angled on side

OM 33-0544 sharp/sharp



DEAVER SCISSORS

5 1/2" (140mm)



OM 34-0142

Straight

OM 34-0140 sharp/blunt
OM 34-0142 sharp/sharp

Curved

OM 34-0150 sharp/blunt
OM 34-0152 sharp/sharp

Most scissors can be serrated upon request.





Straight
OM 53-0330 6" (152mm)
OM 53-0332 7" (179mm)

Curved
OM 53-0331 6" (152mm)
OM 53-0333 7" (179mm)

REYNOLDS SCISSORS

blunt points




OM 34-0336

Straight
OM 34-0334

Curved
OM 34-0336

DOYEN ABDOMINAL SCISSORS

7" (178mm)



OM 34-0258


Straight
OM 34-0256 5 1/2" (140mm)
OM 34-0260 6 3/4" (171mm)
OM 34-0292 9" (229mm)

Curved
OM 34-0258 5 1/2" (140mm)
OM 34-0262 6 3/4" (171mm)
OM 34-0294 9" (229mm)

Serrated
OM 34-0280 6 3/4" (171mm)
OM 34-0282 6 3/4" (171mm)

MAYO SCISSORS

Beveled blades



OM 67-0168




OM 67-0174

Straight
OM 67-0164 sharp/sharp
OM 67-0166 sharp/blunt
OM 67-0168 blunt/blunt

Curved
OM 67-0170 sharp/sharp
OM 67-0172 sharp/blunt
OM 67-0174 blunt/blunt

SIMS SCISSORS

8" (203mm)




OM 34-0327

Straight
OM 34-0325

Curved
OM 34-0327

MAYO-NOBLE SCISSORS

6 1/2" (165mm) Beveled blades



OM 34-0312

Straight
OM 34-0306 5 1/2" (140mm)
OM 34-0310 6 3/4" (171mm)

Curved
OM 34-0308 5 1/2" (140mm)
OM 34-0312 6 3/4" (171mm)

MAYO-STILLE SCISSORS

Rounded blades



OM 34-0338

FERGUSON SCISSORS

7" (178mm)
Angled on flat




OM 34-0354

Straight
OM 34-0348

Curved
OM 34-0354

SISTRUNK SCISSORS

5" (127mm) Delicate




Straight
OM 34-0356

Curved
OM 34-0358

SISTRUNK SCISSORS

5 1/2" (140mm) Heavy pattern



OM 34-0462

Straight


OM 34-0436 4 1/2" (114mm)
OM 34-0458 6" (152mm)
OM 34-0466 7" (178mm)
OM 34-0472 8" (203mm)
OM 34-0476 9" (229mm)
OM 34-0486 11" (279mm)

Curved

OM 34-0438 4 1/2" (114mm)
OM 34-0492 5 3/4" (146mm) delicate
OM 34-0462 6" (152mm)
OM 34-0470 7" (178mm)
OM 34-0494 7" (178mm) delicate
OM 34-0474 8" (203mm)
OM 34-0478 9" (229mm)
OM 34-0488 11" (279mm)

METZENBAUM SCISSORS


blunt points



Angled sideways

OM 34-0502 5 1/2" (140mm)
OM 34-0504 8 1/2" (216mm)

BRAUN EPISIOTOMY SCISSORS



OM 34-0560


OM 34-0562

Curved

OM 34-0560
5" (127mm)
flat blunt tips

OM 34-0562
7" (178mm)
flat blunt tips

RAGNELL SCISSORS



OM 34-0602

OM 34-0607

Straight

OM 34-0602 smooth sharp/sharp

Curved


OM 34-0607 smooth sharp/sharp

Angled

OM 34-0612 smooth sharp/sharp

KELLY FISTULA SCISSORS

6 1/4" (159mm)



OM 53-0072

OM 53-0138

Straight

OM 53-0064 3 1/2" (89mm)

OM 53-0068 4 1/2" (114mm)

OM 53-0080 4 1/8" (105mm) delicate

Curved

OM 53-0066 3 1/2" (89mm)

OM 53-0072 4 1/2" (114mm)


OM 53-0084 4 1/8" (105mm) delicate

Angled on side

OM 53-0138 4 1/2" (114mm)

IRIS SCISSORS

Sharp points



OM 53-0160

OM 53-0164

Straight


OM 53-0160 blunt points

Curved

OM 53-0164 blunt points

STRABISMUS SCISSORS

4" (102mm)



OM 53-0168

OM 53-0174

Straight

OM 53-0168 blunt/blunt

OM 53-0174 sharp/sharp


Curved

OM 53-0172 blunt/blunt

OM 53-0176 sharp/sharp

STEVENS TENOTOMY SCISSORS

4 1/4" (108mm)




Straight

OM 53-0265 3 3/4" (95mm)

Curved

OM 53-0267 4" (102mm)

IRIS SCISSORS RIBBON TYPE




OM 53-0316
3 3/4" (95mm)
slightly curved
sharp/sharp

GRADLE SCISSORS



OM 53-0318
4" (102mm)
slightly curved
ribbon type
sharp/sharp

GRADLE SCISSORS



OM 35-0515
4" (102mm)
angled on flat
sharp outside edges


**BIRO DERMAL
NAEVUS SCISSORS**



OM 35-0078
4 5/8" (117.5mm)
with suture hole


OM 35-0078SC
4 5/8" (117.5mm)
with suture hole
super cut

LITTLER SCISSORS




OM 69-0106
3 1/2" (89mm)
curved stainless

NAIL SCISSORS



OM 35-0080
3 1/2" (89mm)
curved on side
delicate


**SHORTBENT STITCH
SCISSORS**



Straight
OM 35-0068
3 1/2" (89mm)
delicate


OM 35-0092
4 3/4" (121mm)
delicate

**SPENCER STITCH
SCISSORS**




OM 35-0084
4 1/4" (108mm)

**HILL STITCH
SCISSORS**



OM 35-0104
4 3/4" (121mm)
curved on side


NORTHBENT STITCH SCISSORS



OM 35-0120
5 1/2" (140mm)

LITTAUER STITCH SCISSORS

Note: Tungsten carbide (TC) is one of the hardest alloys made today. Our TC scissors have tungsten carbide firmly bonded to both cutting edges. We guarantee that they will stay sharp for at least three years. The finger rings are gold plated for easy identification.



**TC OPERATING
SCISSORS**



OM 35-0600



OM 35-0610

Straight
OM 35-0600
5 1/2" (140mm)
sharp/blunt

Curved
OM 35-0610
5 1/2" (140mm)
sharp/blunt



TC MAYO SCISSORS

Beveled blades



OM 35-0620



OM 35-0622


Straight
OM 35-0620 5 1/2" (140mm)
OM 35-0624 6 3/4" (171mm)
OM 35-0628 9" (229mm)

Curved
OM 35-0622 5 1/2" (140mm)
OM 35-0626 6 3/4" (171mm)
OM 35-0630 9" (229mm)



TC METZENBAUM SCISSORS


Blunt points



OM 35-0646

Straight
OM 35-0640 6" (152mm)
OM 35-0644 7" (178mm)
OM 35-0652 9" (229mm)

Curved
OM 35-0642 6" (152mm)
OM 35-0646 7" (178mm)
OM 35-0654 9" (229mm)




OM 35-0670 OM 35-0672

Straight
OM 35-0670
sharp points

Curved
OM 35-0672
sharp points

TC IRIS SCISSORS

4 1/2" (114mm)




OM 35-0674 OM 35-0676

Straight
OM 35-0674
blunt points

Curved
OM 35-0676
blunt points


TC STRABISMUS SCISSORS

4 1/2" (114mm)



Curved
OM 35-0510
4 1/2" (114mm)
flat tips,
one serrated and
one carbide blade

BLEPHAROPLASTY SCISSORS




OM 31-0692

OM 31-0668 3 1/2" (89mm)
OM 31-0680 4 1/2" (114mm)
OM 31-0692 5 1/2" (140mm)
OM 31-0704 7 1/4" (184mm)
OM 31-0708 8" (203mm) one large finger ring


LISTER BANDAGE SCISSORS

The delicate offset shanks provide better leverage. The fine knob slides easily under finger and toe bandages.



OM 31-0792
5 1/2" (140mm)
shanks angled on side
one serrated blade


HI-LEVEL BANDAGE SCISSORS



For finger and toe bandages.

Straight
OM 31-0804
5 1/2" (140mm)
one serrated blade

KNOWLES BANDAGE SCISSORS



Useful for cutting tape and bandages.

OM 46-0060
8" (203mm)

ESMARCH SHEARS




For heavier plaster.

OM 46-0090
9 1/4" (235mm)
one serrated blade

BRUNS SHEARS

For light plaster and cloth.


With autoclavable plastic finger rings.



OM 31-0830
7 1/2" (191mm)
one serrated blade

OM 31-0832
6" (152mm)
Lister-type
(not pictured)

UTILITY SCISSORS



OM 31-0834
7 1/2" (191mm)
one serrated blade
heavy duty

HERCULES SCISSORS



*Effective for cloth
or cotton wadding.*

OM 31-0836
8" (203mm)
with wire cutting notch
one serrated blade

ORTHOMED UTILITY SHEARS



OM MX 5-700
Pink handle

OM MX 5-702
Blue handle

NURSES SCISSORS

5 1/2" (140mm)

Notch for cutting cerclage wire up to 1mm.

OM 35-0198



OM 35-0198
one serrated blade
angular

OM 35-0200
one serrated blade
angular carbide

WIRE CUTTING SCISSORS

4 3/4" (121mm)




OM 35-0194

Straight
OM 35-0168
smooth

OM 35-0170
one serrated blade

Curved
OM 35-0192
smooth

OM 35-0194
one serrated blade



WIRE CUTTING SCISSORS

4" (102mm)




OM 36-0063



OM 36-0067

Straight
OM 36-0063

Curved
OM 36-0067



HARTMANN MOSQUITO FORCEPS

3 1/2" (89mm)



OM 36-0164



OM 36-0165


Straight
OM 36-0164

Curved
OM 36-0165



MICRO MOSQUITO FORCEPS

5" (127mm) Extremely fine tips



Straight

OM 49-0602	1 1/2" (38mm)
OM 49-0604	2" (51mm)
OM 49-0606	2 1/2" (64mm)
OM 49-0608	3" (76mm)
OM 49-0610	3 1/2" (89mm)

Curved


OM 49-0618	1 1/2" (38mm)
OM 49-0620	2" (51mm)
OM 49-0622	2 1/2" (64mm)
OM 49-0624	3" (76mm)
OM 49-0626	3 1/2" (89mm)

JOHNS HOPKINS BULLDOG CLAMPS



HALSTED MOSQUITO FORCEPS

5" (127mm)





Straight

OM 36-0156
OM 36-0160 delicate
OM 36-0076 1x2 teeth

Curved

OM 36-0158
OM 36-0162 delicate
OM 36-0078 1x2 teeth



Straight



OM 36-0166

Curved

OM 36-0168

KELLY FORCEPS

5 1/2" (140mm)



Straight



OM 36-0196

Curved

OM 36-0198

RANKIN-KELLY FORCEPS

6 1/4" (159mm)



Straight


OM 36-0106

Curved

OM 36-0108


PROVIDENCE HOSPITAL FORCEPS

5 1/2" (140mm)




CRILE FORCEPS

OM 36-0176	OM 36-0178
Straight	Curved
OM 36-0176	OM 36-0178
5 1/2" (140mm)	5 1/2" (140mm)
OM 36-0180	OM 36-0182
5 1/2" (140mm)	5 1/2" (140mm)
extra delicate	extra delicate
OM 36-0186	OM 36-0188
6 1/4" (159mm)	6 1/4" (159mm)



ROCHESTER PEAN FORCEPS


OM 36-0306	OM 36-0308
Straight	
OM 36-0306	5 1/2" (140mm)
OM 36-0316	6 1/4" (159mm)
OM 36-0326	7 1/4" (184mm)
OM 36-0336	8" (203mm)
OM 36-0346	9" (229mm)
OM 36-0356	10" (254mm)
OM 36-0360	12" (305mm)
Curved	
OM 36-0308	1/2" (140mm)
OM 36-0318	6 1/4" (159mm)
OM 36-0328	7 1/4" (184mm)
OM 36-0338	8" (203mm)
OM 36-0348	9" (229mm)
OM 36-0358	10" (254mm)
OM 36-0362	12" (305mm)



ROCHESTER-OGHSNER FORCEPS

1 x 2 Teeth

OM 36-0416	OM 36-0418
Straight	Curved
OM 36-0406	OM 36-0408
5 1/2" (140mm)	5 1/2" (140mm)
OM 36-0416	OM 36-0418
6 1/4" (159mm)	6 1/4" (159mm)
OM 36-0426	OM 36-0428
7 1/4" (184mm)	7 1/4" (184mm)
OM 36-0436	OM 36-0438
8" (203mm)	8" (203mm)
OM 36-0446	OM 36-0448
9" (229mm)	9" (229mm)
OM 36-0456	OM 36-0458
10" (254mm)	10" (254mm)



MIXTER FORCEPS


Fully curved

OM 36-0557
OM 36-0566
OM 36-0557
5 1/4" (133mm)
delicate
OM 36-0562
6 1/4" (159mm)
OM 36-0564
7 1/4" (184mm)
OM 36-0566
9" (229mm)



OM 36-0570
7" (171mm)
serrated

GEMINI MIXTER FORCEPS



Straight
OM 36-0744

Curved
OM 36-0746

ADSON HEMOSTATIC FORCEPS

7 1/4" (184mm)



OM 36-0518


Straight
OM 36-0516
6 1/4" (159mm)

OM 36-0536
8" (203mm)

Curved
OM 36-0518
6 1/4" (159mm)

OM 36-0538
8" (203mm)

ROCHESTER-CARMALT FORCEPS




Straight
OM 36-0816
6 1/2" (165mm)

OM 36-0826
7 1/2" (191mm)

Curved
OM 36-0818
6 1/2" (165mm)

OM 36-0828
7 1/2" (191mm)

FERGUSON ANGIOTRIBE FORCEPS




OM 59-0064 OM 59-0066

Half curved
OM 59-0064

Fully curved
OM 59-0066


SCHNIDT FORCEPS

7 1/2" (191mm)



Straight
OM 36-0766
5 3/4"
(146mm)

PENNINGTON FORCEPS



OM 43-0105

OM 43-0105 4 3/4" (121mm)
4x5 teeth delicate

OM 43-0107 5 3/4" (146mm)
4x5 teeth delicate

OM 43-0110 6" (152mm)
3x4 teeth standard


OM 43-0112 6" (152mm)
4x5 teeth standard

OM 43-0116 6" (152mm)
5x6 teeth standard

OM 49-0628 10"
5x6 teeth standard


OM 49-0629 12"
5x6 teeth standard

ALLIS TISSUE FORCEPS



OM 43-0124
7 1/2" (191mm)
5x6 teeth

ALLIS TISSUE FORCEPS



OM 43-0168
6" (152mm)
3x4 teeth

ALLIS-JUDD TISSUE FORCEPS



OM 43-0310

OM 43-0310
6 1/4" (159mm)

OM 43-0324
8" (203mm)


OM 43-0326
9 1/2" (241mm)

BABCOCK FORCEPS



OM 36-0696
6 1/4" (159mm)
3x3 prongs

LAHEY VULSELLUM FORCEPS




OM 43-0524

Straight
OM 43-0520

Curved
OM 43-0524

DOYEN CLAMP


9" (229mm) Longitudinal serrations



OM 37-0444

OM 37-0444
7" (178mm)
smooth

**PRESBYTERIAN TUBE
OCCLUDING FORCEPS**



OM 37-0448

OM 37-0448
5 1/2" (140mm)

OM 37-0450
7 1/4" (184mm)

OM 37-0452
8" (203mm)

VORSE TUBE OCCLUDING FORCEPS

Serrated




OM 37-0453

OM 37-0453 5 1/2" (140mm)
OM 37-0454 7 1/4" (184mm)
OM 37-0455 8" (203mm)

VORSE TUBE OCCLUDING FORCEPS

Serrated, with guard




Straight
OM 39-0644
8" (203mm)

OM 39-0646
11" (279mm)


Curved
OM 39-0648
8" (203mm)

OM 39-0650
11" (279mm)

**UTILITY AND
STERILIZER FORCEPS**



MAGILL ENDOTRACHEAL CATHETER INTRODUCING FORCEPS



OM 37-0470
2 1/2" (64mm)

OM 37-0478
3 1/4" (83mm)


JONES TOWEL FORCEPS



OM 37-0458
3 1/2" (89mm)

OM 37-0462
5 1/4" (133mm)

BACKHAUS TOWEL FORCEPS



OM 37-0464
5 1/4" (133mm)
with ball stops

BACKHAUS-ROEDER TOWEL CLAMP




OM 47-0466

OM 37-0466
4" (102mm)
non-perforating


OM 37-0468
5 1/4" (133mm)
non-perforating

**LORNA (EDNA)
TOWEL CLAMP**



OM 37-0456
5 3/4" (146mm)
non-perforating

PEERS TOWEL CLAMP



OM 37-0388


Straight
OM 37-0366 7" (178mm)
OM 37-0386 9 1/2" (241mm)

Straight, smooth
OM 37-0370 7" (178mm)
OM 37-0390 9 1/2" (241mm)

Curved, serrated
OM 37-0368 7" (178mm)
OM 37-0388 9 1/2" (241mm)

Curved, smooth
OM 37-0372 7" (178mm)
OM 37-0392 9 1/2" (241mm)


**FOERSTER SPONGE
FORCEPS**



OM 47-0466

OM 38-0076 4 1/2" (114mm)
OM 38-0080 5" (127mm)
OM 38-0084 5 1/2" (140mm)
OM 38-0088 6" (152mm)
OM 38-0094 8" (203mm)
OM 38-0100 10" (254mm)
OM 38-0106 12" (305mm)




DRESSING FORCEPS
Serrated



1x2 teeth
OM 38-0264 4 1/2" (114mm)
OM 38-0268 5" (127mm)
OM 38-0272 5 1/2" (140mm)
OM 38-0276 6" (152mm)
OM 38-0280 8" (203mm)
OM 38-0286 10" (254mm)
OM 38-0290 12" (305mm)

2x3 teeth
OM 38-0292 5" (127mm)
OM 38-0296 5 1/2" (140mm)
OM 38-0299 7" (178mm)
3x4 teeth
OM 38-0304 5 1/2" (140mm)

TISSUE FORCEPS



1x2 teeth
OM 38-0322

2x3 teeth
OM 38-0324

FERRIS-SMITH TISSUE FORCEPS
7" (178mm) with platform



OM 39-0076



OM 39-0080

Serrated

OM 39-0076 5" (127mm)


OM 39-0086 6" (152mm)

1x2 teeth

OM 39-0080 5" (127mm)

OM 39-0088 6" (152mm)

SEMKEN FORCEPS



OM 39-0058



OM 39-0060



OM 39-0062



OM 39-0072

OM 39-0058 micro tip


OM 39-0060 delicate tip

OM 39-0062 standard tip

OM 39-0072 carbide tip

ADSON DRESSING FORCEPS

4 3/4" (121mm) serrated



OM 39-0056



OM 39-0064



OM 39-0066



OM 39-0068



OM 39-0070

1x2 teeth

OM 39-0056 micro tip

OM 39-0063 serrated tip

OM 39-0064 delicate tip

OM 39-0066 standard tip

OM 39-0070 delicate, with tying platform


OM 39-0071 carbide tip

2x3 teeth

OM 39-0068

ADSON TISSUE FORCEPS

4 3/4" (121mm)



OM 39-0057



OM 39-0059



OM 39-0050



OM 39-0052

OM 39-0057 micro 1x2 teeth

OM 39-0059 micro serrated

OM 39-0050 1x2 teeth

OM 39-0052 serrated


ADSON DRESSING AND TISSUE FORCEPS

6" (152mm)



OM 39-0092
4 3/4" (121mm)
9x9 teeth
side grasping

ADSON BROWN TISSUE FORCEPS



OM 39-0098
7" (178mm) 3mm
2x3 teeth

OM 39-0099
7" (178mm) 3mm
2x3 teeth
with carbide inserts

ADSON TISSUE FORCEPS



OM 39-0055
TC 1X2 teeth SERR. 6"

ADSON TISSUE FORCEP



OM 39-0136

OM 39-0128 4 1/2" (114mm)
OM 39-0132 5" (127mm)
OM 39-0136 5 1/2" (140mm)
OM 39-0140 6" (152mm)

DRESSING FORCEPS


Serrated with fluted handles




OM 39-0310 4 1/2" (114mm)
OM 39-0314 5" (127mm)
OM 39-0318 5 1/2" (140mm)
OM 39-0322 6" (152mm)

TISSUE FORCEPS

1x2 teeth with fluted handles



Curved
OM 38-0082
5" (127mm)
fine point serrated



OM 39-0192

OM 39-0192 6" (152mm)
OM 39-0196 8" (203mm)
OM 39-0197 9" (229mm)
OM 39-0198 10" (254mm)



OM 53-0372



OM 53-0374

OM 53-0368 4" (102mm)
straight

OM 53-0372 4" (102mm)
half curved

OM 53-0374 4" (102mm)
fully curved

DRESSING FORCEPS

RUSSIAN TISSUE FORCEPS

EYE DRESSING FORCEPS

4" (102mm) serrated



OM 49-0168



OM 49-0170



OM 49-0172


OM 49-0168 6 1/4" (159mm) 1.5mm
OM 49-0170 8" (203mm) 1.5mm
OM 49-0172 9 1/4" (235mm) 1.5mm
OM 49-0180 6 1/4" (159mm) 2.0mm
OM 49-0182 8" (203mm) 2.0mm
OM 49-0184 9 1/4" (235mm) 2.0mm
OM 49-0210 6 1/4" (159mm) 3.5mm
OM 49-0212 8" (203mm) 3.5mm
OM 49-0216 9 1/4" (235mm) 3.5mm

ADSON TISSUE FORCEPS


4 3/4" (121mm)



OM 53-0404



OM 53-0408



OM 53-0420




OM 53-0422

OM 53-0404 straight serrated
OM 53-0408 curved serrated
OM 53-0420 straight 1x2 teeth
OM 53-0422 curved 1x2 teeth

ADSON DRESSING AND TISSUE FORCEPS

6" (152mm)




OM 53-0382 OM 53-0384

1x2 teeth
OM 53-0380 straight
OM 53-0382 half curved
OM 53-0384 full curved

IRIS FORCEPS

4" (102mm)



Extra narrow
OM 38-0166

DRESSING FORCEPS


5 1/2" (140mm)



1x2 teeth delicate
OM 38-0180 straight
OM 38-0181 curved
OM 38-0182 straight TC

GERALD DRESSING FORCEPS


6 3/4" (171mm)



OM 38-0186 1x2 teeth
OM 38-0184 serrated

CUSHING BAYONET TISSUE FORCEPS

7 1/4" (184mm)



OM 39-0074
8 1/4" (210mm)
straight with serrated tips

ADSON BAYONET FORCEPS

8 1/4" (210mm)



OM 39-0450
#3
4 3/4" (121mm)
narrow fine



OM 39-0451
#3C
4 3/8" (111mm)
narrow fine



OM 39-0452
#4
4 3/8" (111mm)
fine



OM 39-0453
#5
4 3/8" (111mm)
extra fine




OM 39-0455
#6
4 1/2" (114mm)
angled pickup



OM 39-0456
#7
4 1/2" (114mm)
curved fine

SWISS JEWELERS STYLE FORCEPS



OM 39-0386

OM 39-0386
3 1/4" (83mm)

OM 39-0390
4 1/2" (114mm)



OM 39-0366

OM 39-0364
3 1/2" (89mm)


OM 39-0366
4 1/2" (114mm)



OM 39-0410

Straight
OM 39-0410

Curved
OM 39-0412



OM 39-0412

FEILCHENFELD FORCEPS


Serrated

PLAIN SPLINTER FORCEPS

Serrated

CARMALT SPLINTER FORCEPS

4 1/4" (108mm) Serrated




OM 39-0478
6" (152mm)
serrated



OM 39-0534
5" (127mm)
trocar in reversible handle

OM 39-0536
trocar point only



OM 39-0509

OM 39-0511

OM 39-0509
5 1/2" (140mm)
straight

OM 39-0511
5 1/4" (133mm)
45 degree angled

OM 39-0513
5 1/2" (140mm)
90 degree angled

VIRTUS SPLINTER FORCEPS

VANGUARD SPLINTER LIBERATOR

WILLIAMS SPLINTER FORCEPS

Serrated

Note: Our carbide needle holders have TC inserts firmly bonded to the inside of both jaws. We guarantee them for three years. The surface profiles of the inserts have been precision milled to provide the surgeon with a firm, non-slip grip on the needle. Because the carbide is harder than the steel needles, carbide needle holders last much longer than ordinary ones. The finger rings are gold plated for easy identification.



OM 41-0064



OM 41-0066

OM 41-0064
smooth

OM 41-0065
smooth carbide

OM 41-0066
serrated

OM 41-0067
serrated carbide

HALSEY NEEDLE HOLDER

5" (127mm) Serrated



OM 41-0164

OM 41-0166

OM 41-0162
5" (127mm)
serrated

OM 41-0164
5 3/4" (146mm)
serrated

OM 41-0166
5 1/2" (140mm)
serrated carbide

HEGAR BAUMGARTNER NEEDLE HOLDER




DERF NEEDLE HOLDER

4 3/4" (121mm)



BROWN NEEDLE HOLDER

Fine serrations




**WEBSTER
NEEDLE HOLDER**

5" (127mm) Smooth




**COLLIER
NEEDLE HOLDER**




**NIEVERT
NEEDLE HOLDER**

5" (127mm) Smooth, one offset finger ring




MAYO HEGAR NEEDLE HOLDER

Serrated




CRILE-WOOD (DEBAKEY) NEEDLE HOLDER

Serrated




OLSEN-HEGAR NEEDLE HOLDER

Combination needle holder and suture scissors



GILLIES NEEDLE HOLDER

Combination needle holder and suture scissors



OM 42-0148

OM 42-0148
7 1/2" (191mm)


OM 42-0150
10 1/2" (267mm)

Carbide
OM 42-0152
7 1/2" (191mm)

OM 42-0154
10 1/2" (267mm)

SAROT NEEDLE HOLDER

Serrated



OM 42-0257

OM 41-0257
5" (127mm)

OM 41-0264
6" (152mm)

OM 41-0265
7" (178mm)


OM 41-0266
8" (203mm)

OM 41-0267
9" (229mm)

OM 41-0268
10 1/4" (260mm)

RYDER NEEDLE HOLDER

Serrated carbide



OM 42-0257

OM 41-0240
5" (127mm)

OM 41-0242
6" (152mm)


OM 41-0244
7" (178mm)

OM 41-0246
8" (203mm)

OM 41-0248
9" (229mm)

RYDER MICRO NEEDLE HOLDER

Serrated carbide




OM 42-0198

OM 41-0198
5 1/2" (140mm)

Carbide
OM 41-0200
5 1/2" (140mm)

MATHIEU NEEDLE HOLDER

Serrated



OM 42-0186

OM 42-0186
straight

OM 42-0188
curved

Carbide
OM 42-0190
straight

OM 42-0192
curved

CASTROVIEJO MICRO NEEDLE HOLDER

5 1/2" (140mm) smooth, with lock



OM 42-0196


OM 42-0196
straight

OM 42-0198
curved

CASTROVIEJO MICRO NEEDLE HOLDER


5 1/2" (140mm) extra delicate, smooth, with lock

OAL = Overall Length
WL = Working Length
FR = French Gauge (1 french = .33mm)



	OAL	FR	WL
OM 42-0832	5 1/2" (140mm)	3	75mm
OM 42-0834	5 1/2" (140mm)	5	75mm
OM 42-0836	5 1/2" (140mm)	7	75mm


BARON SUCTION TUBE



	OAL	FR	WL
OM 42-0823	7" (178mm)	6	90mm
OM 42-0824	7" (178mm)	7	90mm
OM 42-0825	7" (178mm)	8	90mm
OM 42-0826	7" (178mm)	10	90mm
OM 42-0827	7" (178mm)	11	90mm
OM 42-0830	7" (178mm)	12	90mm
OM 42-0823L	9" (229mm)	6	130mm
OM 42-0824L	9" (229mm)	7	130mm
OM 42-0825L	9" (229mm)	8	130mm
OM 42-0826L	9" (229mm)	10	130mm
OM 42-0827L	9" (229mm)	11	130mm
OM 42-0830L	9" (229mm)	12	130mm

FRAZIER SUCTION TUBE


Angled 30°



	OAL	FR	WL
OM 42-0863	7 1/2" (191mm)	7	85mm
OM 42-0864	7 1/2" (191mm)	9	85mm
OM 42-0866	7 1/2" (191mm)	11	85mm
OM 42-0873	9 1/2" (241mm)	7	130mm
OM 42-0874	9 1/2" (241mm)	9	130mm
OM 42-0876	9 1/2" (241mm)	11	130mm

FRAZIER SUCTION TUBE

Angled 90°



	OAL	FR	WL
OM 42-0810	7" (178mm)	7	90mm
OM 42-0811	7" (178mm)	8	90mm
OM 42-0812	7" (178mm)	9	90mm
OM 42-0813	7" (178mm)	10	90mm
OM 42-0814	7" (178mm)	11	90mm
OM 42-0815	7" (178mm)	12	90mm

FRAZIER SUCTION TUBE


Angled 75°



Tear drop opening allows surgeon to regulate force of suction.

OM 42-0880	6 1/4" (159mm)	4 french
OM 42-0881	6 1/4" (159mm)	5 french
OM 42-0882	6 1/4" (159mm)	6 french
OM 42-0883	6 1/4" (159mm)	7 french
OM 42-0884	6 1/4" (159mm)	8 french
OM 42-0885	6 1/4" (159mm)	9 french
OM 42-0890	7" (178mm)	4 french
OM 42-0891	7" (178mm)	5 french
OM 42-0892	7" (178mm)	6 french
OM 42-0893	7" (178mm)	7 french
OM 42-0894	7" (178mm)	8 french
OM 42-0895	7" (178mm)	9 french
OM 42-0901	8" (203mm)	4 french
OM 42-0902	8" (203mm)	5 french
OM 42-0903	8" (203mm)	6 french
OM 42-0904	8" (203mm)	7 french
OM 42-0905	8" (203mm)	8 french
OM 42-0906	8" (203mm)	9 french

FUKUSHIMA SUCTION TUBE



OM 58-0228
11 1/2" (292mm)
chrome

OM 58-0230
front tip only

YANKAUER SUCTION TUBE



1 french= .33mm

OM 42-0844
23 french curved

OM 42-0852
30 french straight

POOLE SUCTION TUBE

9" (229mm) Chrome



OM 29-0070
for blades 61 through 67

**MINI SCALPEL HANDLE
(BEAVER) #3K 4"**



OM 29-0073

**SCALPEL HANDLE MINI
(BEAVER) #3K 6"**




OM 29-0071

**SCALPEL HANDLE MINI
(BEAVER) #3K 5"**



OM 29-0058
for blades
10 through 15

SCALPEL HANDLE #3



OM 29-0059
mm and cm graduations
for blades 10 through 15

SCALPEL HANDLE #3



OM 29-0060
for blades
20 through 25

SCALPEL HANDLE #4



OM 29-0061
extra long
8 1/2" (216mm)
for blades 10 through 15

SCALPEL HANDLE #3L



OM 29-0062
for blades 10 through 15

SCALPEL HANDLE #7



OM 29-0064
extra fine
for blades 10 through 15


SCALPEL HANDLE #9



OM 42-0608	4 1/2" (114mm)
OM 42-0610	5" (127mm)
OM 42-0612	5 1/2" (140mm)
OM 42-0614	6" (152mm)
OM 42-0616	8" (203mm)
OM 42-0618	10" (254mm)
OM 42-0620	12" (305mm)

PROBE WITH EYE

Malleable stainless



With tongue tie

OM 42-0646	4 1/2" (114mm)
OM 42-0647	5" (127mm)
OM 42-0648	5 1/2" (140mm)
OM 42-0649	6" (152mm)

With probe and tongue tie

OM 42-0654	5" (127mm)
OM 42-0655	5 1/2" (140mm)
OM 42-0656	6" (152mm)


GROOVED DIRECTORS



OM 28-0460
7" (178mm) small

OM 28-0468
8" (203mm) large

TAYLOR PERCUSSION HAMMER



OM 28-0492
7 3/4" (197mm)
chrome with brush and needle

BUCK NEUROLOGICAL HAMMER

UTILITY CLIP MINI

OM CK-10

UTILITY CLIP LARGE

OM CK-40



REVERSE SYNOVECTOMY TISSUE FORCEPS

OM 38-0325


7" 1x2 teeth

OM 38-0326

7" 2x3 teeth

The new Synovectomy Forcep clasps onto dense soft tissue. Its reverse design functions like a clothespin, gripping tissue securely. The self-retaining design mitigates any thumb pain surgeons may experience with traditional forceps.





OM 42-0558 blunt left
OM 42-0560 blunt right
OM 42-0562 sharp left
OM 42-0564 sharp right

DESCHAMPS LIGATURE CARRIER

8"(203mm)



OM 53-0450
0.12mm
1x2 teeth

OM 53-0456
0.5mm
1x2 teeth

CASTROVIEJO MICRO TYING FORCEPS

4"(102mm) Straight with tying platform

*Bipolar forceps
Available in most styles and tip sizes.*



OM 68-0026
4"(102mm)
8mm platform
0.5mm tip

JEWELLERS FORCEP