

MALLETS, TAMPS AND IMPACTORS

SECTION K


CONTENTS

Mallets	1-6
Tamps	7-8
Impactors and Drivers	8-10



Mallets exert force on osteotomes, chisels and gouges.
Tamps wedge bone graft into place.

Note: Lead filled mallets should only be used to strike flat surfaces. Do not use to hit sharp objects as doing so can cause lead to flake off.



OM 49-0765
7" (178mm)
solid stainless steel
weight: 4 oz. (120 grams)
head weight: 3 oz. (80 grams)
head diameter: 20mm


MINI Mallet



OM 49-0761
7 1/2" (191mm)
solid stainless steel
weight: 7 oz. (200 grams)
head weight: 3 oz. (80 grams)
head diameter: 19mm


OM 46-0660
8" (203mm)
brass head
stainless steel handle
weight: 9 oz. (260 grams)
head weight: 5 oz. (150 grams)
head diameter: 25mm

LUCAE BONE Mallet



OM 46-0656
7 3/4" (195mm)
lead filled and faced
weight: 7 oz. (200 grams)
head weight: 5 oz. (150 grams)
head diameter: 22mm

Mallet #49



OM 46-0661
7 1/2" (191mm)
lead filled and faced
weight: 10 oz. (285 grams)
head weight: 7 oz. (200 grams)
head diameter: 25mm

GERZOG BONE Mallet



Useful for hand surgery and in areas where access is limited.

Autoclavable to 290° F (143° C).

OM 49-0764
8" (203mm)
with black plastic handle
weight: 7 oz. (200 grams)
head weight: 3 oz. (80 grams)
head diameter: 18mm x 18mm
(square contact surface)
7mm (round contact surface)

NARROW TIP Mallet



OM 49-0760
8" (203mm)
solid stainless steel
weight: 7 oz. (200 grams)
head weight: 5 oz. (150 grams)
head diameter: 30mm


COLLIN Mallet



Autoclavable to 250° F (121° C).

OM 46-0662
7 1/2" (191mm)
with two replaceable nylon caps
anodized aluminum handle
weight: 7 oz. (200 grams)
head weight: 4 oz. (120 grams)
head diameter: 25mm


Mallet
Lightweight



Autoclavable to 250° F (121° C).

OM 46-0657
7 1/2" (191mm)
with replaceable nylon caps
weight: 11 oz. (315 grams)
head weight: 6 oz. (170 grams)
head diameter: 20mm

MEAD Mallet



Autoclavable to
250° F (121° C).

OM 46-0686
7 1/2" (191mm)
weight: 9 oz. (260 grams)
head weight: 5 oz. (150 grams)
head diameter: 26mm

OM 46-0687
7 1/2" (191mm)
weight: 12 oz. (340 grams)
head weight: 8 oz. (240 grams)
head diameter: 30mm

COMBINATION Mallet

With one nylon cap and one solid stainless side



OM 55-0400
7 1/2" (191mm)
weight: 8 oz. (240 grams)
head weight: 4 oz. (120 grams)
head diameter: 19mm
ball peen diameter: 16mm

OM 55-0401
9" (229mm)
weight: 1 lb. (450 grams)
head weight: 8 oz. (240 grams)
head diameter: 25mm
ball peen diameter: 22mm

CLOWARD-STYLE Mallet

Solid stainless steel



OM 46-0656
7 3/4" (195mm)
lead filled and faced
weight: 7 oz. (200 grams)
head weight: 5 oz. (150 grams)
head diameter: 22mm

ORTHOPEDIC Mallet

Short handle, lightweight



*Head is brass filled
and should only be
used to strike flat
surfaces.*

OM 46-0658
7 3/4" (195mm)
weight: 1 lb. 7 oz. (650 grams)
head weight: 1 lb. 4 oz. (570 grams)
head diameter: 38mm

KIRK Mallet



OM 49-0762
7 1/2" (191mm)
solid stainless steel

weight:
2 lbs. (900 grams)

head weight:
1 lb. 5 oz. (600 grams)
head diameter: 38mm

ORTHOPEDIC Mallet

Heavyweight, short handle



OM 49-0766
9" (229mm)
aluminum alloy handle

weight:
15 oz. (420 grams)

head weight:
9 oz. (260 grams)
head diameter: 30mm

BERGMAN Mallet



OM 49-0767
9 1/2" (241mm)
aluminum alloy handle

weight:
1 lb. 1 oz. (470 grams)

head weight:
11 oz. (310 grams)
head diameter: 45mm

BERGMAN Mallet



OM 46-0663
7 1/2" (191mm)
aluminum handle
weight: 2 lbs. (900 grams)
head weight: 1 lb. 4 oz. (570 grams)
head diameter: 40mm

HEATH Mallet

Heavyweight



*Also known as dead
blow mallet.*

*Grooved phenolic handle
for comfortable grip.
Autoclavable to 250° F (121° C).*

OM 46-0690
11 1/2" (292mm)
weight: 1 lb. 5 oz. (600 grams)
head weight: 10 oz. (285 grams)
head diameter: 32mm
dead blow mallet with black plastic
handle, one nylon cap and one
stainless steel side

REPROCUSSION FREE Mallet



OM 46-0667
10 1/2" (267mm)
solid stainless steel
weight: 3 lbs. (1360 grams)
head weight: 2 lbs. (900 grams)
head diameter: 48mm

ORTHOPEDIC Mallet

Heavyweight



*Grooved phenolic handle
for comfortable and safe grip.
Autoclavable to 290° F (143° C).*

OM 46-0688
10" (254mm)
weight: 1 lb. 11 oz. (765 grams)
head weight: 1 lb. (455 grams)
head diameter: 35mm

BLACK PLASTIC HANDLE Mallet



OM 46-0666
10 3/4" (273mm)
solid stainless steel
weight: 2 lbs. (900 grams)
head weight: 1 lb. (455 grams)
head diameter: 32mm

ORTHOPEDIC Mallet

Lightweight



Grooved phenolic handle for comfortable and safe grip. Autoclavable to 290°F (143°C).

OM 46-0689
11" (280mm)
weight: 2 lbs. 9 oz. (1160 grams)
head weight: 2 lbs. (900 grams)
head diameter: 45mm

BLACK PLASTIC Mallet



3 lb Mallet

OM 46-0697
11" (280mm)
9" (229mm)
knurled handle


5 lb Mallet

OM 46-0698
11" (280mm)
9" (229mm)
knurled handle



OM 49-0770
9 3/4" (248mm)

WAGNER Mallet



OM 49-0771
1 pair of metal disks



OM 49-0772
1 pair of plastic disks



OM 49-0773
wrench




OM 12500
9" (230mm)

BONE Mallet


Weight: 1lb. 7oz (700 gr)
Head diameter: 35mm




BONE CHIP PACKER
10" (254mm)




OM 49-0751
2mm x 9mm



OM 49-0752
6mm x 9mm



OM 49-0750
3mm x 16mm



OM 49-0753
6mm x 19mm



OM 49-0754
1/4" diameter



BONE TAMP
Endoscopic





OM 55-0676
17" (432mm)
stainless steel handle
3/8" diameter
.16" tip length
12.59" shaft length
4 1/4" handle length
.280" (7.1mm) shaft diameter



OM 46-0893	2mm diameter
OM 46-0894	3mm diameter
OM 46-0895	4mm diameter
OM 46-0896	5mm diameter
OM 46-0897	6mm diameter
OM 46-0898	8mm diameter
OM 46-0899	10mm diameter

BONE TAMP

With cross serrated tip
6" (152mm)




OM 69-0980P
replaceable nylon cap

Designed to tap bone graft or bone parts into place with minimal trauma to the bone. The compactor head is nylon with a grid surface to prevent slippage.

Autoclavable to 250° (121°C).

OM 69-0980
6 1/2" (165mm)
with replaceable nylon caps
anodized aluminum handle
head diameter: 20mm

BONE IMPACTOR



OM 55-0370	8 1/2"	Impactor, standard offset
OM 55-0371	8 1/2"	Impactor, central canal
OM 55-0372	8 1/2"	Impactor, lateral narrow
OM 55-0373	7"	Impactor, micro offset
OM 55-0374	7"	Impactor, micro lateral
OM 55-0375	8 1/2"	Graft Driver, lumbar, guarded
OM 55-0376	8 1/2"	Graft Driver, anterior, guarded
OM 55-0377	8 1/2"	Graft Driver, lumbar, unguarded
OM 55-0378	8 1/2"	Graft Driver, anterior, unguarded

LUMBAR AND CERVICAL IMPACTORS AND DRIVERS

Serrated tips



OM 55-0384	Bone Graft Holder / Impactor 6" (152mm)	Serr 3.2mm 1/8"
OM 55-0385	Bone Graft Holder / Impactor 6" (152mm)	Serr 4.8mm 3/16"
OM 55-0386	Bone Graft Holder / Impactor 6" (152mm)	Serr 6.4mm 1/4"
OM 55-0387	Bone Graft Holder / Impactor 6" (152mm)	Serr 7.9mm 5/16"
OM 55-0388	Bone Graft Holder / Impactor 10" (254mm)	Serr 3.2mm 1/8"
OM 55-0389	Bone Graft Holder / Impactor 10" (254mm)	Serr 4.8mm 3/16"
OM 55-0390	Bone Graft Holder / Impactor 10" (254mm)	Serr 6.4mm 1/4"
OM 55-0391	Bone Graft Holder / Impactor 10" (254mm)	Serr 7.9mm 5/16"

BONE TAMP

With cross serrated tip
6" (152mm)