


MEASURING DEVICES

SECTION W

CONTENTS

Screw Depth Gauge	1
Calipers	1-2
Rulers	3
Easy Clean Depth Gauge	4



To determine dimensions during surgery.




TOWNLEY CALIPER

OM 48-0314

4" (102mm)

takes inside and outside measurements up to 4" (10cm)



BONE COMPASS

OM 48-0302

8" (203mm)

Used to determine screw length required during internal fracture fixation (mm graduated). Gauges snap apart for easy cleaning.



MINI SCREW DEPTH GAUGE

OM 47-0934

4 3/4" (121mm)

up to 30mm



SMALL SCREW DEPTH GAUGE

OM 47-0936

6 1/4" (159mm)

up to 50mm



LARGE SCREW DEPTH GAUGE

OM 47-0940

8 3/4" (222mm)

up to 100mm




X-LARGE SCREW DEPTH GAUGE

OM 47-0935


11" (279mm)

up to 120mm




CASTROVIEJO CALIPER

OM 53-0040
3 1/2" (89mm)
graduated to 20mm
for hand surgery




CASTROVIEJO CALIPER

OM 53-0042
6 1/2" (165mm)
graduated to 40mm




POLK FINGER GONIOMETER

OM 48-0310
6" (152mm)
cm and inch graduated



ALUMINUM CALIPER

OM 48-0312
9 1/8" (232mm)
inch and mm graduated
up to 6" (180mm)
inside and outside measuring
faces and depth bar



NEURO-CALIPER

OM 48-0315
9" (229mm)
stainless steel
inch and mm graduated
up to 5 3/4" (145mm)
inside and outside measuring
faces and depth bar



OM 39-0692 6" (152mm)
OM 39-0694 8" (203mm)
OM 39-0698 12" (305mm)
OM 39-0700 20" (508mm)

RULER


flexible
inch and mm graduated



Ruler for x-rays.


X-RAY SCALE

OM 48-0299
50cm length
mm graduated



MEASURING TAPE

OM 48-0313
linen, inch and
cm graduated
up to 150cm



Front Back

Reverse side has handy gauge to
measure K-Wires, Steinmann Pins,
twist drills, and pins and rods.

Note:

Gauge Measures	
K-Wires	.028"-.062" .7mm-1.6mm
Pins	5/64"-3/16" 2.0mm-4.76mm
Rods	.197"-1/4" 5.0mm-6.35mm

K-Wire RULER & PIN GAUGE

OM 39-0710
6" (152mm)
inch and mm graduated



EASY CLEAN DEPTH GAUGE

	Screw Size	Hook Probe	Measuring Range
OM 24-0116	1.2-1.5mm	0.8mm	0-30mm
OM 24-0118	2.0-2.4mm	1.2mm	0-30mm
OM 24-0120	2.0-2.4mm	1.2mm	0-50mm
OM 24-0124	2.7mm	1.7mm	0-50mm
OM 24-0128	3.5-4.0mm	2.2mm	0-50mm
OM 24-0130	3.5-4.0mm	2.2mm	0-80mm