

MISCELLANEOUS INSTRUMENTS

SECTION Z

CONTENTS

Lead Hands	1
Forceps Holders	2
Finger Ring Cutter	2
Pin Wheel	2



Note: Oxidation of lead surfaces is a normal process.

**LEAD HAND**

OM 49-0092
9 1/2" (241mm)
child

**LEAD HAND**


OM 49-0090
14" (356mm)
adult

To prolong the life of lead hands, autoclave temperatures should not exceed 248°F (120°C). Avoid using chemical sterilization as they can cause severe corrosion.

Thinner and more flexible model. Longer fingers wrap on ends eliminating side tabs.

**LEAD HAND**

OM 49-0093
9 1/2" (241mm)
adult
2mm thick



ORTHOMED FORCEPS HOLDER

OM 39-0672
6" (152mm)

OM 39-0674
12" (305mm)

OM 39-0673
8" (204mm)




BUNT FORCEPS HOLDER

OM 39-0670



FINGER RING CUTTER

OM 29-0430
1 1/2" (165mm) chrome
with saw blade



WARTENBERG NEUROLOGICAL PIN WHEEL

OM 28-0514 7" (178mm)