

TRAYS AND RACKS

SECTION AA



To hold instruments for sterilization and storage.

CONTENTS

Hibbs-Brun Rack	1
Lambotte Rack	1
Smith Peterson Rack	1
Brun Rack	2
Hoke Rack	2
Cobb Rack	2
OrthoMed K-Wire and Steinmann Pin Sterilization Case	3-6



HIBBS-BRUN RACK

OM 98-0180

12" x 8" x 3" (305mm x 203mm x 76mm)

Holds twelve 9" (229mm) OrthoMed Brun Curettes, Hibbs Chisels, Gouges or Osteotomes



LAMBOTTE RACK

OM 98-0190

11 1/2" x 8" x 3" (292mm x 203mm x 76mm)

Holds twelve 9" (229mm) OrthoMed Lambotte Osteotomes



SMITH PETERSON RACK

OM 98-0195

10" x 10" x 1 7/8"

(254mm x 254mm x 48mm)

Holds seven 8" (203mm) OrthoMed Smith Peterson Osteotomes or Gouges



BRUN RACK

OM 98-0160
 9 1/2" x 8" x 3"
 (241mm x 203mm x 76mm)
 Holds 12 instruments



HOKE RACK

OM 98-0150
 8" x 5 1/2" x 1 1/4"
 (203mm x 140mm x 32mm)
 Holds six instruments



COBB RACK

OM 98-0170
 13" x 12" x 1 7/8"
 (330mm x 305mm x 48mm)
 Holds nine 11" OrthoMed Cobb currettes,
 elevators, or gouges



Designed to provide easy access and identification.

ORTHOMED K-Wire AND STEINMANN PIN STERILIZATION CASE

OM 71-0010

11 3/4" x 5 1/4" x 2" (closed)

(298mm x 133mm x 51mm)

The OrthoMed K-Wire and Steinmann Pin Sterilization Case holds popular 4" (102mm), 6" (152mm) and 9" (229mm) K-Wires and 9" (229mm) Steinmann Pins.

Accommodates four each of the following K-Wire diameters:

.028" (.7mm)
.035" (.9mm)
.045" (1.1mm)
.062" (1.6mm)

Accommodates two each of the following Steinmann Pin diameters:

5/64" (2.0mm)
3/32" (2.4mm)
7/64" (2.8mm)
1/8" (3.2mm)
9/64" (3.5mm)
5/32" (4.0mm)
.177" (4.5mm)

- Removable lid fits under case – taking up less room in the O.R.

- Nylon coated brackets prevent metal-to-metal contact – protecting sharp ends.

- Contents are held in place when closed – preventing shifting.

- Secure locking side brackets help prevent accidental opening.


- Handles at each end for easy transport.

- Anodized aluminum case is lightweight and built to withstand abuse.

- Fully autoclavable to 270°F (132°C).



Case in closed position.



ORTHOMED K-Wire AND STEINMANN PIN STERILIZATION CASE

OM 71-0009

13 5/8" x 5 1/4" x 2 1/4" (closed)

(346mm x 133mm x 57mm)

The OrthoMed K-Wire and Steinmann Pin Sterilization Case holds popular 4" (102mm), 6" (152mm) and 9" (229mm) K-Wires and 9" (229mm) Steinmann Pins.

Accommodates four each of the following K-Wire diameters:

.028"	(.7mm)
.035"	(.9mm)
.045"	(1.1mm)
.062"	(1.6mm)

Accommodates two each of the following Steinmann Pin diameters:

5/64"	(2.0mm)
3/32"	(2.4mm)
7/64"	(2.8mm)
1/8"	(3.2mm)
9/64"	(3.5mm)
5/32"	(4.0mm)
.177"	(4.5mm)

- Removable lid fits under case – taking up less room in the O.R.
- Nylon coated brackets prevent metal-to-metal contact – protecting sharp ends.
- Contents are held in place when closed – preventing shifting.
- Secure locking side brackets help prevent accidental opening.
- Handles at each end for easy transport.
- Anodized aluminum case is lightweight and built to withstand abuse.
- Fully autoclavable to 270°F (132°C).



Case in closed position.

The OrthoMed K-Wire sterilization container holds 4" (102mm) or 6" (152mm) K-Wires.

Accommodates six each of the following K-Wire diameters:

.028"	(.7mm)
.035"	(.9mm)
.045"	(1.1mm)
.054	(1.4mm)
.062"	(1.6mm)

- K-Wire sizes are clearly noted on the inside of case for quick identification and reordering.

- Lid folds back and acts as a support to hold the case in an upright position – taking up less room in the O.R.

- Top of case locks when closed helping to prevent accidental opening.

- Case has a one piece design – eliminates the possibility of losing or switching case parts.

- Anodized aluminum case is lightweight and built to withstand abuse.

- Fully autoclavable to 270°F.



Case in open upright position.

ORTHOMED K-Wire STERILIZATION CONTAINER

OM 71-0008

6 1/4" x 4 1/8" x 3/4" (closed)

(158mm x 105mm x 19mm)



Case in closed position.



ORTHOMED K-Wire AND STEINMANN PIN STERILIZATION CASE

OM 71-0011

13 5/8" x 5 1/4" x 2 1/4" (closed)

(346mm x 133mm x 57mm)

The OrthoMed, Inc. K-Wire and Steinmann Pin Sterilization Case holds popular 9" K-Wires and Steinmann Pins.

Accommodates twelve each of the following K-Wire diameters:

.035"	(.9mm)
.045"	(1.1mm)
.062"	(1.6mm)

Accommodates four each of the following Steinmann Pin diameters:

5/64"	(2.0mm)
3/32"	(2.4mm)

Accommodates two each of the following Steinmann Pin diameters:

7/64"	(2.8mm)
1/8"	(3.2mm)
9/64"	(3.5mm)
5/32"	(4.0mm)

- Removable lid fits under case – taking up less room in the O.R.
- Nylon coated brackets prevent metal-to-metal contact – protecting sharp ends.
- Contents are held in place when closed – preventing shifting.
- Secure locking side brackets help prevent accidental opening.
- Handles at each end for easy transport.
- Anodized aluminum case is lightweight and built to withstand abuse.
- Fully autoclavable to 270°F (132°C).



Case in closed position.